

GINGKO

G50 Motor-Driven Tissue Grinder



G50 electric tissue grinder is a brushless motor drive small, convenient, all-digital control handheld tissue grinder / homogenizer, in a small centrifuge tube grinding tiny protein, DNA and tissue particles play grinding and mixing role.

Features

3000-8000 RPM adjustable speed

Built-in constant speed control by external load remains constant speed output 4-bit LCD screen, buttons to set the speed Using the sensorless brushless gun technology do not need to rely on any speed sensor to calculate the real-time speed.

The loop current control technology, to prevent the load operation of the motor burned effectively protect the battery and power, High torque brushless motor as a power output, the output power is the same size as similar products 4-5 times, Brushless motor greatly extend the life of the motor.

Technical parameters

Power supply: 12V DC 5A (vehicle power supply or battery box).

Speed range: 3000 ~ 8000rpm.

Motor Control: Brushless motor with constant speed control system.

Output Power: 50W continuous output (instantaneous output 70W).

Dimensions: 170 × 50 × 26mm

Product Weight: 145g

