

# ALL-AMERICAN® PORTABLE STERILIZERS

## Electric Models

Wisconsin Aluminum Foundry was first in the field to develop an immersion-type pressure sterilizer for physicians' offices, hospitals and laboratories. The electric ALL-AMERICAN® Sterilizer is a complete, self-contained unit with an immersion heating element, automatic thermostatic control, automatic release valve, and pilot light. It connects to any outlet, is ready for use instantly, and ensures complete, efficient sterilization.

- Available in 120V or 240V
- Operational on 50Hz or 60Hz
- Calrod immersion heating element
- 3-wire grounded cord and plug for safety
- Optional support bases elevate electric sterilizers above table or counter surface to prevent heat damage
- UL approved
- One year warranty
- Units include aluminum inner container with rack
- Interior stainless steel support stand
- Made in U.S.A.

### SPECIFICATIONS

#### Model 75X (41 qt/39 liter)

**Model 75X-120:** 120 Volt, 50/60 Hz 1650 watts/13.75 amps

**Model 75X-240:** 240 Volt, 50/60 Hz 1650 watts/6.88 amps

**Gross Capacity** 41 qt / 39 liter

Overall Height 19" / 48.3cm

Bottom Height 14 1/4" / 36.2cm

Inside Diameter 15 1/4" / 38.7cm

Unit Weight 45 lbs. / 20.4kg.

#### Inner Container No. 4156

Height 10 1/4" / 26cm

Diameter 14" / 35.6cm

Circumference 44 1/2" / 113cm

Capacity 27.3 qt / 25.8 liter

Volume 1578in<sup>3</sup> / 25,856cm<sup>3</sup>

Carton Dimensions 24" x 24" x 21"

61cm x 61cm x 53.3cm

Shipping Weight 51 lb. / 23.1kg.

Unit Pack: 1 Cube: 7



#### 75X Features:

- 64% greater capacity compared to 25X
- 120 or 240 volt operation
- 1650 watt heating element for faster operation
- More accurate TPI thermocouple, snapswitch and pilot light
- 60 minute mechanical bell timer
- Extra large 41 qt/39 liter capacity