

Viscosight Bath 30L



FEATURES

1. Microprocessor PID Control for 0.1°C Precision Temperature Control from Ambient +5°C to 99°C.
2. Front mount glass window for sample view.
3. Timer, alarm, Auto-tuning function.
4. Built in safety and convenience
5. Powerful circulation pumps with 30liter/min, capacity enhances temperature uniformity.
6. Adjustable stainless steel shelf

Specifications:

Model		LVB-630SP	LVB-650SP
Model		LVB-630S	LVB-650S
Controller		Digital PID controller	
Dimension (wxdxh) mm	Inner	320x300x300	500x300x350
	Outer	573x400x380	725x400x430
Capacity		30 Liter	52 Liter
Temperature	Range	Ambient +5°C to 99°C	
	Set Accuracy	±0.1°C	
	Uniformity	±0.2° C	
Heater		1.5kW	2kW
Controller		Digital PID Multi function controller	
Display		LED 4 Digital Display	
Timer		99hr 59 min/Continuous	
Pump		Magnetic drive Max Flow: 30 Liter/min Max Height: 4.3 Meter	
Material	Interior	Stainless Steel (SUS304)	
	Exterior	Powder coated steel	
	Window	Vacuum sealed tempered safety glass	
Safety		Over current and leakage breaker , Over temp , Cut-Off	
Electric Supply		110V, 60HZ	