

Midi plus Horizontal Electrophoresis System, ME7-7-10

Product Introduction

Our Midi series have been designed for routine horizontal gel electrophoresis with wide degrees of versatility. All of our gel tanks are made from injection molded construction for durable and leak proof environment. Varieties of special combs are available to maximize your sample size. Midi plus includes the 70x70 and 70x100mm gel trays for a maximum of 64 samples, designed for low to medium numbers of samples.

FEATURES

- Low buffer consumption
- Single molded tank
- Two tray options available (ME7-7-10, ME10-7-10)
- Easy sample loading
- Multi-channel pipette compatibility
- Great indications for gel making and running
- No tape, clamps or springs needed



| | |
|------------------------------|-------------------------------|
| Cat. No. | ME7-7-10 |
| Dimension(mm)(WxLxH) | 210x90x90 |
| Gel dimension(mm)(WxL) | 70x70 70x100 |
| Maximum sample each tray(mm) | 32 for 70x70 64 for 70x100 |
| Rapid casting gel | Use casting dams |

SPECIFICATION:

| | |
|-------------------------|--|
| Cat. No. | ME7-7-10 |
| Dimension (WxLxH) | 8.3"x3.5"x3.5" (210x90x90mm) |
| Gel Dimension (WxL) | 2.8"x2.8" (70x70mm) 2.8"x3.9" (70x100mm) |
| Maximum Sample | 32 for 2.8"x2.8" (70x70mm) Tray (2x16 sample combs) |
| | 64 for 2.8"x3.9" (70x100mm)Tray (4x16 sample combs) |
| Buffer Volume | 0.06gal (225ml) |
| Construction | Injection molded construction |
| | Durable, leak proof environment for complete safety and long life |
| Cassette Type Electrode | Inexpensive, easy to replace |
| | Made of 99.99% corrosion resistant pure platinum |
| Electrical Safety | Lid can only be fastened in one way |
| | Upon lid removal, power is disconnected from buffer chamber |
| Rapid Casting Gel | Use casting dams |